



T:1-800-265-9830 F:1-877-636-8944

1000 Series Bed

TROUBLESHOOTING Guide

If you experience trouble with your bed, try the following solutions.

1. Problem: The Back Rest up & down operation does not work

- 1.1. Verify if the power cord is plugged into the wall outlet properly.
 - 1.1.1. If answer is yes, go to step 1.2.
 - 1.1.2. If answer is no, plug power cord into wall outlet, then go to step 1.2.
- 1.2. Verify if the power indicating light on control box is on.
 - 1.2.1. If answer is yes, go to step 1.3.
 - 1.2.2. If answer is no, hold the reset button for 10 seconds until the light turns on, then go to step 1.3.
- 1.3. Verify if the back rest motor power cord is plugged into the control box.
 - 1.3.1. If answer is yes, go to step 1.4.
 - 1.3.2. If answer is no, plug back rest motor power cord into the control box, then go to step 1.4.
- 1.4. Verify if the back rest operation is functional/operates by remote control (Do not lie on the bed)
 - 1.4.1. If answer is yes, the back rest up & down operation should function normally.
 - 1.4.2. If answer is no, please contact customer service.

2. Problem: The Foot Rest up & down operation does not work

- 2.1. Verify if the power cord is plugged into the wall outlet properly.
 - 2.1.1. If answer is yes, go to step 2.2.
 - 2.1.2. If answer is no, plug the power cord into wall outlet. Then go to step 2.2.
- 2.2. Verify if the power indicating light on control box is on.
 - 2.2.1. If answer is yes, go to step 2.3.
 - 2.2.2. If answer is no, hold the reset button for 10 seconds until the light turns on, then go to step 2.3.
- 2.3. Verify if the foot rest motor power cord is plugged into the control box.
 - 2.3.1. If answer is yes, go to step 2.4.
 - 2.3.2. If answer is no, plug foot rest motor power cord into the control box, then go to step 2.4.
- 2.4. Verify if the foot rest is functional/operates by remote control (Do not lie on the bed)
 - 2.4.1. If answer is yes, the foot rest up & down operation should function normally.
 - 2.4.2. If answer is no, please contact customer service.

3. Problem: The Bed horizontal up & down operation does not work

- 3.1. Verify if objects or articles are obstructing the up & down function of the bed.
 - 3.1.1. If answer is no, go to step 3.2.
 - 3.1.2. If answer is yes, remove the obstructing objects and or articles then go to step 3.2.
- 3.2. Verify if the power cord is plugged into the wall outlet properly.
 - 3.2.1. If answer is yes, go to step 3.3.
 - 3.2.2. If answer is no, plug the power cord into the wall outlet, then go to step 3.3.
- 3.3. Verify if the power indicating light on the control box is on.
 - 3.3.1. If answer is yes, go to step 3.4.
 - 3.3.2. If answer is no, hold the reset button for 10 seconds until the light turns on, then go to step 3.4.
- 3.4. Verify if the lifting motor power cord is plugged into the control box.
 - 3.4.1. If answer is yes, go to step 3.5.
 - 3.4.2. If answer is no, plug the lifting motor power cord into the control box, then go to step 3.5.
- 3.5. Verify if the up & down mobility operation is functional/operates by remote control (Do not lie on the bed)
 - 3.5.1. If answer is yes, the up & down mobility operation should function normally.
 - 3.5.2. If answer is no, please contact customer service.

4. Problem: Bed Trendelenburg operation does not work

- 4.1. Verify if objects or articles are obstructing the up & down function of the bed.
 - 4.1.1. If answer is yes, remove the obstructing objects and or articles then go to step 4.2.
 - 4.1.2. If answer is no, go to step 4.2.
- 4.2. Verify if the power cord is plugged into the wall outlet properly.
 - 4.2.1 If answer is yes, go to step 4.3.
 - 4.2.2 If answer is no, plug power cord into wall outlet, then go to step 4.3.
- 4.3. Verify if the power indicating light on control box is on.
 - 4.3.1. If answer is yes, go to step 4.4.
 - 4.3.2. If answer is no, hold the reset button for 10 seconds until the light turns on, then go to step 4.4.
- 4.4. Verify if the lifting motor power cord is plugged into the control box.
 - 4.4.1. If answer is yes, go to step 4.5.
 - 4.4.2. If answer is no, plug the lifting motor power cord into the control box, then go to step 4.5.
- 4.5. Verify that the Trendelenburg operation is functional/operates by remote control (Do not lie on the bed)
 - 4.5.1. If answer is yes, your Trendelenburg operation should function normally.
 - 4.5.2. If answer is no, please contact customer service.

5. Problem: Side rail up & down operation does not work

- 5.1. Verify if objects or articles are obstructing the movement of the side rails.
 - 5.1.1. If answer is no, go to step 5.2.
 - 5.1.2. If answer is yes, remove the obstructing objects and or articles then go to step 5.2.
- 5.2. Verify if locking devices fully open and close.
 - 5.2.1. If answer is yes, your side rail up & down operation should function normally.
 - 5.2.2. If answer is no, please contact customer service.

6. Problem: Bed will not roll

- 6.1. Verify around and underneath the bed for any objects or articles that may be obstructing the movement of the bed.
 - 6.1.1. If answer is no, go to step 6.2.
 - 6.1.2. If answer is yes, remove the obstructing objects and or articles then go to step 6.2.
- 6.2. Verify if the casters are on brake position.
 - 6.2.1. If answer is yes, release brake.
 - 6.2.2. If answer is no, please contact customer service.

If the troubleshooting guide did not solve your problems, please contact one of our company's customer service representatives for further assistance:



T: 1-800-265-9830 F: 1-877-636-8944

www.advancedhealthcare.ca

2221- 46th Ave., Lachine, Quebec H8T 3C9