EASyS

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Dear EASyS user,

The EASyS rehabilitation stroller that you purchased is an innovative product made by the company Thomashilfen. After several years' of research we have developed a new concept for rehabilitation strollers' that is one of the safest and most advanced seating systems available.

In order to make the operation of the EASyS seating system and frame safe, practical, and comfortable for you and especially for your child, please read this user manual completely before use.

Should you have any further questions or issues, please contact your specialist supplier or Thomashilfen directly.

For more information about Thomashilfen and our products, please visit our website at www.thomashilfen.com

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Contents

| | Page |
|------------------------------------|------|
| General safety instructions..... | 12 |
| Technical data | 13 |
| Setting up and assembly | 13 |
| Adjustment options..... | 13 |
| Mounting the belts | 14 |
| Care instructions..... | 14 |
| Terms of guarantee / CE mark | 15 |

Please observe the following safety instructions:

- We recommend that you adjust the EASyS seat growth settings every 3 months to best fit your child. If necessary, let your therapist / prosthetist / dealer help you.
- This stroller is intended to transport only one child. More than one child in or on the stroller may result in injury.
- Check all control elements and fastening screws monthly. The sliding bars should be kept free of sand.
- Please remember to re-tighten the screws and clamping levers (hip angle adjustment) after each re-setting.
- Always lock the stroller brakes before placing or removing a child from the stroller; and before adding any accessories or performing maintenance.
- **CAUTION:** wet or flat wheels may affect the braking effect. If not in use or stopping for a longer period, engage the stroller foot brake in order to prevent accidental rolling.
- Do not exceed the maximum size and weight capacity of this stroller. Please refer to the literature for specification details.
- Leaning or reaching far too any side can cause the stroller to be unstable and tip over. Always use the seat belt in addition to close observation to prevent this type of movement.
- Attaching bags, equipment, or other heavy items may cause the stroller to become unstable. These items are best stowed in the basket / equipment tray below the seat unit.
- We recommend that you always secure your child with the 2-point pelvic belt (standard equipment) or another retention belt (available as an accessory). Never leave your child unattended in the stroller; even if the brakes are engaged and the seat belt is in use. Be especially careful about other children playing around the stroller. The stroller should only be operated by a mature individual.
- This stroller contains many small pieces that can be a choking hazard. Never leave children unattended in or around the stroller. Choking hazards include, but are not limited to: packaging materials, screws, nuts, and plastic caps which can become loose or dislodged with use. Inspect the stroller frequently to insure safety and proper operation.
- The pelvic belt (seat belt) supplied with the stroller, as standard equipment, is NOT intended for use as a safety belt for transportation in a motor vehicle. When transporting the child seated in the stroller in a motor vehicle, always ensure that the tie down and occupant restraint system meet all current regulations.
- When using an anterior trunk support (H, butterfly, or vest style harness) always use a pelvic belt and/or an abduction block to prevent sliding that can result in strangulation.
- Always be alert when pushing a child in the stroller. Give particular attention to using the stroller when on or near slopes, uneven ground, stairs, obstacles and when near moving motor vehicles. Stairs should only be attempted with two or more capable attendants.
- Strolling after dark can be extremely hazardous. Please use great caution, obey all traffic laws, use lights and/or light clothing to increase visibility to others.
- When working with the seat unit, please observe the recommendations in the user manual, especially when dismantling it and folding it together.
- Any modifications that alter the frame, such as drilling holes in the frame, will void the approved use of this stroller for transportation in a vehicle.
- Before use, always check that the seat is securely fastened to the frame and that the recline and tilt mechanisms are securely fastened.

Tip: More information is available on our website www.thomashilfen.com

Technical data

| EASyS seating system | EASyS 1 | EASyS 2 |
|--|-------------------------------------|--|
| Seat depth | 19 - 30 cm / 7.5 – 11.8" | 27 - 38 cm / 10.6 – 15" |
| Seating width | 19 - 31 cm / 7.5 – 12.2" | 25 - 36 cm / 9.8 – 14.2" |
| Back height | 40 - 65 cm / 15.8 – 25.6" | 53 - 73 cm / 20.9 – 28.8" |
| Lower leg length | 17 - 30 cm / 6.7 – 11.8" | 25 - 36 cm / 9.8 – 14.2" |
| Knee angle (adjustable in steps of 15°) | 90° - 180° | 90° - 180° |
| Hip angle | 90° - 180° | 90° - 180° |
| Foot angle | 78° to 100° | 78° to 100° |
| Seat tilt-in-space on J/A underframe | -15° to +45° | -10° to +40° |
| Size of footrest (w x d) | 29,8 x 17,5 cm / 11.7 x 6.8" | 33,5 x 21,5 cm / 13.5 x 8.5" |
| Overall dimensions (w x h x l) | 38 x 98 x 64 cm / 15 x 38.6 x 25.2" | 45 x 116 x 76 cm / 17.7 x 45.7 x 29.9" |
| Folded dimensions (w x h x l) | 38 x 29 x 88 cm / 15 x 11.4 x 34.7" | 45 x 40 x 81 cm / 17.7 x 15.8 x 31.9" |
| Carrying capacity | 35 kg / 77.1 lb | 40 kg / 88.1 lb |
| Weight | 7.7 kg / 17 lb | 9.4 kg / 20.7 lb |
| Seat turning range | 180° | 180° |

Setting up and assembly

The footrest is to be fixed to the knee joints. Place the grid joints one upon the other, and connect them with the screws and the washers (coming from the inside) and wing nuts (coming from the outside). (fig. 1)

Mounting the lateral upholsteries

Put the lateral backrest upholstery over the lateral tube and velcro it to the back plate (the snap fastener points inwards). Close the Velcro of the lateral backrest upholstery completely. Attach the lateral upholstery to the lateral backrest upholstery with the snap fastener. In addition, velcro it to the back plate. (fig. 2-4)

Seat upholstery

Fix the two lower Velcro fasteners around the tubes of the footrest.

Inserting the seat unit

Place the seat unit on the rear crossbar and engage forwards. The EASyS click-system engages automatically. (fig. 5)

Removal of the seat unit

To remove the seat unit, push the EASyS click-lever under the seat forwards and take the seat unit off the chassis forwards.

Adjustment options

Seat width

Loosen the hip guide screws from the seat top surface (under the seat upholstery) with an Allen key. Move the hip guide to the appropriate width and retighten the screws. (fig. 6)

Seat depth

Loosen the front Allen screws (also for the hip guide setting) on the surface and pull out the seat plate. Retighten the screws when the seat depth is correct. (fig. 7)

Tip: For easy adjustment of the seat depth grip the frontal edge of the seat with one hand and press against the bottom part of the backrest (concentric) with the other hand.

Back height

Loosen both black knobs to adjust the backrest variably to the desired height. (fig. 9)

Lateral back support

The tube bow is infinitely adjustable in height by loosening both knurled nuts at the upper cross tube. The maximum height is achieved when the locking pin clicks in audibly. (fig. 10). Press the locking pin (only on the right-hand side) to reduce the height of the lateral support. (fig. 11)

Tip: The extended lateral support offers additional protection to children who are primarily transported in a lying position.

Backrest adjustment

Loosen the backrest clamping lever and move the backrest to the desired position. The seat unit automatically slides forwards on its rails, shifting its centre of gravity to maintain stability in all seat positions. (fig. 12)

The EASyS clamping lever for hip angle adjustment is equipped with an extra safety locking pin. It must be loosened during adjustment. After adjustment, make sure that the locking pin is engaged again. (fig. 13)

Lower leg length

Release the locking device by loosening the black knobs on the right and the left of the bottom of the footrest; the footrest can be freely adjusted. (fig. 14)

Foot angle

Loosen the 4 screws on the right and left of the rear side of the footrest. Adjust the desired foot angle and firmly re-tighten the screws. (fig. 15)

Knee angle

Set the knee angle by loosening the side wing screws (on both sides). (fig. 16)

Tip: The footrest can be flipped up to facilitate getting on and out of the stroller.

Mounting the belts

Cut the backrest padding in the quilted seam depending on the desired mounting height or open the backrest padding flap. Pull the zip down to open the padding, unhook and release the Velcro fastener.

Pass the two upper belts through the holes in the upholstery and insert them into the front back plate. (fig. 17).

To facilitate mounting pull the back plate out slightly. For this, loosen the two knurled screws at the back (fig. 9)

Insert the belts into the cam buckles and adjust the desired belt length (fig. 18). Reclose the backrest padding flap with the Velcro fastener. Adjust the desired back height and retighten the knurled screws.

The two bottom belt ends can be attached with cam buckles to either the back plate or the plastic element at the side (fig. 19-20). Guide the belts left and right next to the back upholstery to the corresponding cam buckle and fix them.

Tip: The belt length can also be adjusted from the front. Pull the D-ring to shorten the belt.

The front belts of the 4-point pelvic belt are attached right and left at the plastic element with a cam buckle, as well. Guide the belt through the slot in the width plate, insert it into the cam buckle, and fix it (fig. 21).

Care instructions

Upholstery

The upholstery can be easily removed. The Trevira CS material and the interlaced spacing yarns are made of high-quality cloth which can be washed at up to 40°C. The material is printed with non-poisonous dyes which are free of heavy metals. Please be aware that the high-quality materials may also become bleached as a result of long-term and intensive sunlight or due to frequent washing. The upholstery is flame resistant (BS 5852 part 1). (fig. 22)

Velcro fasteners

In order to maintain the functionality of the Velcro fasteners, brush them down from time to time. Always close the Velcro fasteners when washing them.

Terms of guarantee / CE mark

Thomashilfen is certified in accordance with the ISO standard DIN EN ISO 13485:2003 and guarantees the quality of the products they supply. Both material selection and processing meet the very highest of manufacturing standards. Every single part is always subjected to final inspection.

Thomashilfen offers you a guarantee of 3 years for all frame parts, commencing on the date of purchase. This guarantee incorporates all claims which influence functionality. Damages due to improper use (e.g. overloading), and natural wear and tear are excluded.

EASyS meets the requirements of the European standards EN 12182 and 12183 and is provided with a CE marking.

Enjoy using your EASyS rehabilitation stroller!

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