

ThevoChair- The Mobility Chair

Order Form – 2011 Retail

Date	Customer #	PO#	Patient
Customer			
Address			
Contact	Phone ()	Fax ()	



Back with MiS Micro-Stimulation®

Improves body perception of the person with dementia. This intensifies the positive aspect of the rocking motion.

Gentle rocking

Requires little effort thanks to special pendulums to gently move body joints.

Revolving platform

Rotates by 360° offering a range of views and activities and the ease of sitting down and getting up without risk.

Foot stool (option)

Works in harmony with the ThevoChair and provides additional movement impulses.

Safety

Constructed for optimal stability while rocking. The rocking runners move gently, while the chair remains stable.

Fabric

The breathable fabric can be wiped clean and disinfected. It is moisture-resistant and fire-retardant.

Benefits:

- Dementia patients become more receptive and less irritated
- Possibility to reduce or eliminate medication
- Individuals attain better balance
- Rocking movement promotes exercise
- Exceptional design provides added stability and prevents tipping

ThevoChair

Product	Item no.	Retail Price
ThevoChair	<input type="checkbox"/> EX63900	\$ 3,399.00
Cover for Base of Chair	<input type="checkbox"/> EX63912	\$374.00
Foot Stool	<input type="checkbox"/> EX63910	\$624.00
Incontinence Covers for Chair	<input type="checkbox"/> EX63913	\$374.00

Specs

ThevoChair	
Height of back	72 cm / 28.4"
Height (total, with case)	110 cm / 43.3"
Height (without case)	108 cm / 42.5"
Breadth (total)	73 cm / 28.7"
Depth of seat	47 cm / 18.5"
Breadth of seat	55 cm / 18.5"
Back angle	112° / 115°
Seat height (with case)	48 cm / 18.9"
Seat height (without case)	46 cm / 18.1"
Max. load-bearing capacity	140 kg / 308.6 lbs
Foot Stool (accessory)	
Height (max)	41 cm / 16.1"
Breadth	44 cm / 17.3"
Depth	44 cm / 17.3"

This product manufactured by Thomashilfen carries a 36 month limited warranty on all frame parts, commencing on the date of purchase.