

siis
Therapy Chair



Models

Thomashilfen

Therapy Chair SiiS size 1

size 2



size 1

Areas of application:

- Cerebral apraxia
- Cognitive processing deficiencies / ADD
- Tonus dysregulation deficiencies

Item code	8300	8305
Seat height	42 - 50 cm / 16.5 - 19.7"	42 - 61 cm / 16.5 - 24"
Seat width	19 cm / 7.5"	19 cm / 7.5"
Height adjustment of impulse generators	18 cm / 7.1"	18 cm / 7.1"
Iliac crest support		
- height adjustment	8 cm / 3.2"	8 cm / 3.2"
- depth adjustment	3,5 cm / 1.4"	3,5 cm / 1.4"
Total dimensions (w x h x d)	54 x 88 x 54 cm 21.3 x 34.7 x 21.3"	60 x 88 x 56 cm 23.6 x 34.7 x 22"
Push bar height	88 cm / 34.7"	88 cm / 34.7"
Wheel size	10 cm / 3.9"	10 cm / 3.9"
Wheel distance front	53 cm / 20.9"	58 cm / 22.8"
Wheel distance rear	54 cm / 21.3"	55 cm / 21.7"
Max. weight capacity	50 kg / 110 lb	120 kg / 265 lb
Weight	13 kg / 29 lb	14 kg / 31 lb
Body size	max. 135 cm / max. 53.2"	from 130 cm on / from 51.2" on

Therapy Chair SiiS

plus ⊕



Areas of application:

- Cerebral apraxia
- Cognitive processing deficiencies / ADD
- Tonus dysregulation deficiencies
- Perfectly suitable for the use in sheltered workshops

Item code	8306
Seat height	49 - 71 cm / 19.3 - 28"
Seat width	22 cm / 8.7"
Height adjustment of impulse generators	18 cm / 7.1"
Iliac crest support	
- height adjustment	8 cm / 3.2"
- depth adjustment	3,5 cm / 1.4"
Total dimensions (w x h x d)	60 x 88 x 56 cm 23.6 x 34.7 x 22"
Push bar height	88 cm / 34.7"
Wheel size	10 cm / 3.9"
Wheel distance front	58 cm / 22.8"
Wheel distance rear	55 cm / 21.7"
Max. weight capacity	120 kg / 265 lb
Weight	14 kg / 31 lb
Body size	from 140 cm / 55.1" on

Therapy Chair Siis

mini



Areas of application:

- Cerebral apraxia
- Cognitive processing deficiencies / ADD
- Tonus dysregulation deficiencies
- Suitable for the use in the kindergarten

Item code	8301
Seat height	32 - 38 cm / 12.6 - 15"
Seat width	15 cm / 5.9"
Impulse generators	
- height adjustment	27,5 cm / 10.8"
- depth adjustment	3,5 cm / 1.4"
Iliac crest support	
- height adjustment	8 cm / 3.1"
- depth adjustment	3,5 cm / 1.4"
Total dimensions (w x h x d)	54 x 50 x 68 cm 21.3 x 19.7 x 26.8"
Wheel size	10 cm / 3.9"
Wheel distance front	53 cm / 20.9"
Wheel distance rear	54 cm / 21.3"
Max. weight capacity	22 kg / 49 lb
Weight	11 kg / 24 lb
Body size	max. 125 cm / 49.2"

Therapy Chair Siis

high



Areas of application:

- Cerebral apraxia
- Cognitive processing deficiencies / ADD
- Tonus dysregulation deficiencies
- Features: less abduction
- Sheltered workshops

Item code	8307
Seat height	53 - 80 cm / 20.9 - 31.5"
Seat width	19 cm / 7.5"
Impulse generators	
- height adjustment	27,5 cm / 10.8"
- depth adjustment	3,5 cm / 1.4"
Iliac crest support	
- height adjustment	8 cm / 3.1"
- depth adjustment	3,5 cm / 1.4"
Total dimensions (w x h x d)	60 x 88 x 56 cm 23.6 x 34.7 x 22"
Push bar height	88 cm / 34.7"
Wheel size	10 cm / 3.9"
Wheel distance front	58 cm / 22.8"
Wheel distance rear	55 cm / 21.7"
Max. weight capacity	120 kg / 265 lb
Weight	15 kg / 33 lb
Body size	from 145 cm / 57.1" on

Therapy Chair SiiS

narrow



Areas of application:

- Cerebral apraxia
- Cognitive processing deficiencies / ADD
- Tonus dysregulation deficiencies
- Features: less abduction

Item code	8308
Seat height	42 - 61 cm / 16.5 - 24"
Seat width	16 cm / 6.3"
Height adjustment of impulse generators	18 cm / 7.1"
Iliac crest support	
- height adjustment	8 cm / 3.2"
- depth adjustment	3,5 cm / 1.4"
Total dimensions (w x h x d)	60 x 88 x 56 cm 23.6 x 34.7 x 22"
Push bar height	88 cm / 34.7"
Wheel size	10 cm / 3.9"
Wheel distance front	58 cm / 22.8"
Wheel distance rear	55 cm / 21.7"
Max. weight capacity	120 kg / 265 lb
Weight	14 kg / 31 lb
Body size	from 135 cm / 53.2" on