

# ThevoActiv – The mobility bed

## Thevo – The mobility bed


The two components underframe and mattress provide for a soft and at the same time dynamic lying of the patient. They automatically transfer every slightest movement of the patient, even the movement of breathing, into micro size counter movements.

## ThevoActiv

- ✓ the sleep of the patient is not being disturbed
- ✓ if required, actors can be driven in order to give movement impulses to the patient
- ✓ suitable for high-risk patients
- ✓ pressure sore prevention and treatment (up to and including grade III acc. to EPUAP)

### Mattress

Special air channels in the mattress cut provide for a soft lying comfort and an optimum ventilation of the skin.



### Wing suspensions

They take in micromotions and return them to the patient.

### Actors

Provide for pressure relief at endangered skin areas.


### Flexible side element

Fits every hospital bed.

### Control unit

Operates and controls the actors. A simple menu selection offers a variety of movement programmes. This provides for a stimulation individually adapted to the patient.

### Programmes:



## ModelChoice

90 x 200 cm / 35.4 x 78.7"    100 x 200 cm / 39.4 x 78.7"

Version I:			90 x 200 cm / 35.4 x 78.7"	100 x 200 cm / 39.4 x 78.7"
body weight 30 - 50 kg / 66.1 – 110.2 lb	ThevoActiv 50	incl. jersey cover	65102	65101
	ThevoActiv 50	incl. incontinence cover	65112	65111
	ThevoActiv 100	incl. jersey cover	65122	65121
Version II:				
body weight 50 - 100 kg / 110.2 – 220.5 lb	ThevoActiv 100	incl. incontinence cover	65132	65131
	ThevoActiv 135	incl. jersey cover	65142	65141
Version III:				
up to a body weight of 135 kg / 297.6 lb	ThevoActiv 135	incl. incontinence cover	65152	65151

More sizes available on request!

## ThevoActiv cover

ThevoActiv can be used either with a jersey or an incontinence cover.



Jersey cover

Incontinence cover

### Shared features:

- sewn zipper
- washable at 65°C / 149°F
- suitable for tumble dryers
- autoclavable
- suitable for the re-use

### Special features:


- Jersey**
  - good microclimate
- Incontinence cover**
  - breathable
  - cleaning by wiping disinfection

For further information please contact your local dealer.



ThevoActiv

Thevo  
by Thomashilfen



Nursing scientific consultation by



Effect proven and checked by:



Thevo  
by Thomashilfen

Walkmühlenstraße 1  
D - 27432 Bremervörde · Germany  
Phone: +49 (0) 47 61 / 886-63 or -68  
Fax: +49 (0) 47 61 / 886-19  
www.thevo-division.com

Thevo – The mobility bed