



Recovery5 Cover Cleaning Instructions

ROUTINE CARE: Wipe Clean with neutral soapsuds and lukewarm water and allow cover to completely air dry before putting sheets and blankets back on.

DO NOT USE HARSH CLEANERS OR SOLVENTS.

LAUNDERING IS NOT RECOMMENDED

DISINFECTION: Most phenolic or quaternary type disinfectants can be used. User must dilute disinfectants and germicides in accordance with manufacturer's instructions. It is recommended that all cleaning agents first be spot tested on an inconspicuous part of the fabric to ensure that they do not damage the fabric.

Important: DO NOT LEAVE ON SURFACE TO DRY. This will burn the top layer of the fabric.

BLOOD CONTAMINATION: Contamination with blood can be disinfected with a 1:10 dilution of household bleach (5.25% sodium hypochlorite) as recommended by the CDC (US Department of Health and Human Services, February 1989).

IMPORTANT

- Iodophor type disinfectants (Betadine, for example) will stain fabrics.
- Prolonged contact with bleach solutions, improper dilution ratios for cleaning agent, or poor rinsing of fabric surface after disinfecting and cleaning may cause damage to both the appearance and performance of this fabric.
- Do not leave any cleaning agents (i.e. Virox) on surface of Mattress to dry

Recovery5 Cover Cleaning Instructions

ROUTINE CARE: Wipe Clean with neutral soapsuds and lukewarm water and allow cover to completely air dry before putting sheets and blankets back on.

DO NOT USE HARSH CLEANERS OR SOLVENTS.

LAUNDERING IS NOT RECOMMENDED

DISINFECTION: Most phenolic or quaternary type disinfectants can be used. User must dilute disinfectants and germicides in accordance with manufacturer's instructions. It is recommended that all cleaning agents first be spot tested on an inconspicuous part of the fabric to ensure that they do not damage the fabric.

Important: DO NOT LEAVE ON SURFACE TO DRY. This will burn the top layer of the fabric.

BLOOD CONTAMINATION: Contamination with blood can be disinfected with a 1:10 dilution of household bleach (5.25% sodium hypochlorite) as recommended by the CDC (US Department of Health and Human Services, February 1989).

IMPORTANT

- Iodophor type disinfectants (Betadine, for example) will stain fabrics.
- Prolonged contact with bleach solutions, improper dilution ratios for cleaning agent, or poor rinsing of fabric surface after disinfecting and cleaning may cause damage to both the appearance and performance of this fabric.
- Do not leave any cleaning agents (i.e. Virox) on surface of Mattress to dry