



ADVANCED HEALTH CARE PRODUCTS G.P.

Enhancing Your Quality of Life

www.advancedhealthcare.ca

PRODUCT GUIDE



the ROHO group
shape fitting technology™

MATTRESS



Caring for Canadians



the ROHO group
shape fitting technology™

SELECTAIR™

sleepmatterzzz®

a ROHO group product



the ROHO group
shape fitting technology™

CUSHIONS



ADVANCED HEALTH CARE PRODUCTS G.P.

Enhancing Your Quality of Life

www.advancedhealthcare.ca

Phone: 514-636-7575 Toll Free: 800-265-9830 Fax: 877-636-8944

Please refer to the sizing chart for guidance in selecting the proper cushion size.

BACK SYSTEMS

JetStream Pro™ with CF Technology

Extremely light weight back utilizing Carbon Fiber technology. Available in three back heights and five widths.

JETSTREAM PRO STANDARD BACK SYSTEM ADJUSTABLE 16" / 41cm high

JSPI4SFA - JSPI8SFA in 1" increments

JETSTREAM PRO MID BACK SYSTEM ADJUSTABLE 13" / 33cm high

JSPI4MFA - JSPI8MFA in 1" increments

JETSTREAM PRO LOW BACK SYSTEM ADJUSTABLE 10" / 25cm high

JSPI4LFA - JSPI8LFA in 1" increments

CFTFC Replacement Cover (specify height and width)

CFTFC Replacement Foam & Cover (specify size)

JSPPRIVACYFIAP JetStream Pro Privacy Flap



• All systems include back support, mounting hardware, foam pad cushion, cover, installation wrench and instruction manual.

JetStream Pro™ Lateral Support

- | | |
|--|---|
| <input type="checkbox"/> JSPLATSTDSWRTL | JSP STD-LAT Swing Away Right Large Pad |
| <input type="checkbox"/> JSPLATSTDSWLFLG | JSP STD-LAT Swing Away Left Large Pad |
| <input type="checkbox"/> JSPLATSTDSTRTLG | JSP STD-LAT Stationary Right Large Pad |
| <input type="checkbox"/> JSPLATSTDSTLFLG | JSP STD-LAT Stationary Left Large Pad |
| <input type="checkbox"/> JSPLATSTDRTAC | JSP STD-LAT Stationary Right Active Pad |
| <input type="checkbox"/> JSPLATSTDLFAC | JSP STD-LAT Stationary Left Active Pad |
| <input type="checkbox"/> JSPLATLOWSWLG | JSP LOW-LAT Swing Away Large Pad |
| <input type="checkbox"/> JSPLATLOWSTLG | JSP LOW-LAT Stationary Large Pad |
| <input type="checkbox"/> JSPLATLOWSTAC | JSP LOW-LAT Stationary Active Pad |



+ Requires use of the JetStream Pro Adjustable Model

RETROBACK™ Back Support System

Light weight back support provides rigidity, while increasing sitting tolerance and reducing lower back fatigue

- | | |
|-----------------------------------|---------------------------|
| Adjustable Model | Width of Wheelchair |
| <input type="checkbox"/> RETRO16C | 15" to 16" / 38cm to 41cm |
| <input type="checkbox"/> RETRO18C | 17" to 18" / 43cm to 46cm |
| <input type="checkbox"/> RETRO20C | 19" to 20" / 48cm to 51cm |



SEAT CUSHION PRODUCTS

HIGH-PROFILE® Therapeutic Cushions with Covers.

Recommended for individuals with a history of skin breakdown; 4" cell height.

SINGLE COMPARTMENT - STANDARD ADULT SIZES

- IR88C 15"w x 15"d x 4"h (38cm x 38cm x 10cm)
- IR89C 15"w x 17"d x 4"h (38cm x 43cm x 10cm)
- IR99C 17"w x 17"d x 4"h (43cm x 43cm x 10cm)
- IR910C 17"w x 19"d x 4"h (43cm x 47cm x 10cm)
- IR109C 19"w x 17"d x 4"h (47cm x 43cm x 10cm)
- IR1010C 19"w x 19"d x 4"h (47cm x 47cm x 10cm)
- C-ROHO High profile "Custom Cushion" (specify size)
- C-ROHO "Custom Design Cushion" (provide diagram)

PEDIATRIC SIZES

- IR66C 12"w x 12"d x 4"h (30cm x 30cm x 10 cm)
- IR67C 12"w x 13"d x 4"h (30cm x 33cm x 10 cm)



DUAL COMPARTMENT - This cushion is divided into two separate air compartments. Can be used for positioning, amputees, scoliosis, etc.

ADULT SIZES

- 2R88C 15"w x 15"d x 4"h (38cm x 38cm x 10 cm)
- 2R89C 15"w x 17"d x 4"h (38cm x 43cm x 10 cm)
- 2R99C 17"w x 17"d x 4"h (43cm x 43cm x 10 cm)
- 2R910C 17"w x 19"d x 4"h (43cm x 47cm x 10 cm)
- 2R109C 19"w x 17 d x 4"h (47cm x 43cm x 10cm)
- 2R1010C 19"w x 19 d x 4"h (47cm x 47cm x 10cm)
- C-ROHO High profile "Custom" Cushion (specify size)
- C-ROHO "Custom Design" Cushion (provide diagram)

PEDIATRIC SIZES

- 2R66C 12"w x 12"d x 4"h (30cm x 30cm x 10 cm)
- 2R67C 12"w x 13"d x 4"h (30cm x 33cm x 10 cm)

LOW-PROFILE® Therapeutic Cushions with Covers.

Recommended for active individuals without skin problems; 2" cell height.

SINGLE COMPARTMENT - STANDARD

ADULT SIZES

- IR88LPC 15"w x 15"d x 2"h (38cm x 38cm x 5cm)
- IR89LPC 15"w x 17"d x 2"h (38cm x 43cm x 5cm)
- IR99LPC 17"w x 17"d x 2"h (43cm x 43cm x 5cm)
- IR910LPC 17"w x 19"d x 2"h (43cm x 47cm x 5cm)
- IR109LPC 19"w x 17"d x 2"h (47cm x 43cm x 5cm)
- IR1010LPC 19"w x 19"d x 2"h (47cm x 47cm x 5cm)
- C-ROHO Low-Profile "Custom" (specify size)
- C-ROHO "Custom Design" Cushion (provide diagram)

PEDIATRIC SIZES

- IR66LPC 12"w x 12"d x 2"h (30cm x 30cm x 5 cm)
- IR67LPC 12"w x 13"d x 2"h (30cm x 33cm x 5cm)


DUAL COMPARTMENT - This cushion is divided into two separate air compartments. Can be used for positioning, amputees, scoliosis, etc.

ADULT SIZES

- 2R88LPC 15"w x 15"d x 2"h (38cm x 38cm x 5cm)
- 2R89LPC 15"w x 17"d x 2"h (38cm x 43cm x 5cm)
- 2R99LPC 17"w x 17"d x 2"h (43cm x 43cm x 5cm)
- 2R910LPC 17"w x 19"d x 2"h (43cm x 47cm x 5cm)
- 2R109LPC 19"w x 17"d 2"h (47cm x 43cm x 5cm)
- 2R1010LPC 19"w x 19"d x 2"h (47cm x 47cm x 5cm)
- C-ROHO High profile "Custom" Cushion (specify size)
- C-ROHO "Custom Design" Cushion (provide diagram)

PEDIATRIC SIZES

- 2R66LPC 12"w x 12"d x 2"h (30cm x 30cm x 5cm)
- 2R67LPC 12"w x 13"d x 2"h (30cm x 33cm x 5 cm)



ROHO® MID-PROFILE THERAPEUTIC Cushions with Covers.

Recommended for active individuals at increased risk for skin breakdown.

SINGLE COMPARTMENT - STANDARD

ADULT SIZES

- IR77MPC 14" w X 14" d X 3" h (35 cm x 35 cm x 7cm)
- IR87MPC 15" w X 14" d X 3" h (38 cm x 35 cm x 7cm)
- IR88MPC 15" w X 15" d X 3" h (38 cm x 38 cm x 7cm)
- IR89MPC 15" w X 17" d X 3" h (38 cm x 43 cm x 7cm)
- IR99MPC 17" w X 17" d X 3" h (43 cm x 43 cm x 7cm)
- IR910MPC 17" w X 19" d X 3" h (43 cm x 47 cm x 7cm)
- IR109MPC 19" w X 17" d X 3" h (47 cm x 43 cm x 7cm)
- IR108MPC 19" w X 15" d X 3" h (47 cm x 43 cm x 7cm)
- IR1010MPC 19" w X 19" d X 3" h (47 cm x 47 cm x 7cm)
- C-ROHO Mid Profile "Custom" Cushion (please specify suze)

PEDIATRIC SIZES

- IR66MPC 12" w X 12" d X 3" h (30 cm x 30 cm x 7cm)
- IR67MPC 12" w X 13" d X 3" h (30 cm x 33 cm x 7cm)

DUAL COMPARTMENT - STANDARD

ADULT SIZES

- 2R88MPC 15" w X 15" d X 3" h (38 cm x 38 cm x 7cm)
- 2R89MPC 15" w X 17" d X 3" h (38 cm x 43 cm x 7cm)
- 2R99MPC 17" w X 17" d X 3" h (43 cm x 43 cm x 7cm)
- 2R910MPC 17" w X 19" d X 3" h (43 cm x 47 cm x 7cm)

PEDIATRIC SIZES

- 2R66MPC 12" w X 12" d X 3" h (30 cm x 30 cm x 7cm)
- 2R67MPC 12" w X 13" d X 3" h (30 cm x 33 cm x 7cm)



QUADTRO® SELECT™ -HIGH PROFILE Cushions with Covers. FEATURING ISO FLO™ MEMORY CONTROL™ the unique four separate sections of the cushion allows for progressive positioning for short and long term changes.

Standard cell height: 4". Available in 3" and 1" cell height as a custom option.

- QS88C 15" w X 15" d X 4" h (38 cm x 38 cm x 10cm)
- QS89C 15" w X 17" d X 4" h (38 cm x 43 cm x 10cm)
- QS99C 17" w X 17" d X 4" h (43 cm x 43 cm x 10cm)
- QS109C 19" w X 17" d X 4" h (47 cm x 43 cm x 10cm)
- QS910C 17" w X 19" d X 4" h (43 cm x 47 cm x 10cm)
- QS1010C 19" w X 19" d X 4" h (47 cm x 47 cm x 10cm)
- QSC Custom Cushion (please specify size)
- QSC "Custom Design" Cushion (provide diagram)

PEDIATRIC SIZES

- QS66C 12" w X 12" d X 4" h (30 cm x 30 cm x 10cm)
- QS67C 12" w X 13" d X 4" h (30 cm x 33 cm x 10cm)



QUADTRO® SELECT™ -LOW PROFILE Cushions with Covers. FEATURING ISO FLO™ MEMORY CONTROL™ ideal for active user. 2" cell height and four compartments promote lateral, interior and posterior stability.

Standard cell height: 2". Available in 3" and 1" cell height as a custom option.

- QS88LPC 15" w X 15" d X 2" h (38 cm x 38 cm x 5 cm)
- QS89LPC 15" w X 17" d X 2" h (38 cm x 43 cm x 5 cm)
- QS99LPC 17" w X 17" d X 2" h (43 cm x 43 cm x 5 cm)
- QS910LPC 17" w X 19" d X 2" h (43 cm x 47 cm x 5 cm)
- QS109LPC 19" w X 17" d X 2" h (47 cm x 43 cm x 5 cm)
- QS1010LPC 19" w X 19" d X 2" h (47 cm x 47 cm x 5 cm)
- QSC Custom Cushion (please specify size)
- QSC Custom Design Cushion (provide diagram)

PEDIATRIC SIZES

- QS66LPC 12" w X 12" d X 2" h (30 cm x 30 cm x 5 cm)
- QS67LPC 12" w X 13" d X 2" h (30 cm x 33 cm x 5 cm)



QUADTRO® SELECT™ MID PROFILE Cushions with Covers. Featuring ISO FLO™ Memory Control™ that offers shape fitting capabilities while the user is seated, allowing quick and easy on-demand adjustment to maximize function. 3" cell height.

ADULT SIZES

- QS88MPC 15" w x 15" d x 3" h (38 cm x 38 cm x 5 cm)
- QS89MPC 15" w x 17" d x 3" h (38cm x 43cm x 5cm)
- QS99MPC 17" w x 17" d x 3" h (43cm x 43cm x 5cm)
- QS910MPC 17" w x 19" d x 3" h (43cm x 47cm x 5cm)
- QS109MPC 19" w x 17" d x 3" h (47cm x 43cm x 5cm)
- QS1010MPC 19" w x 19" d x 3" h (47cm x 47cm x 5cm)
- QSC Custom Cushion (please specify size and 2" or 4" perimeter)




HIGH-PROFILE® LTC Cushions with Covers.

Recommended for individuals requiring ease of cleaning and heavy usage. Ideal in a Long Term Care Setting. 4" Cell Height.

ADULT SIZES


- IR88LTC 15" w x 15" d x 4" h (38cm x 38cm x 10cm)
- IR89LTC 15" w x 17" d x 4" h (38cm x 43cm x 10cm)
- IR99LTC 17" w x 17" d x 4" h (43cm x 43cm x 10cm)
- IR910LTC 17" w x 19" d x 4" h (43cm x 47cm x 10cm)
- IR109LTC 19" w x 17" d x 4" h (47cm x 43cm x 10cm)
- IR1010LTC 19" w x 19" d x 4" h (47cm x 47cm x 10cm)
- C-ROHO High profile LTC "Custom Cushion" (specify size)
- C-ROHO "Custom Design LTC Cushion" (provide diagram)



LOW-PROFILE® LTC Cushions with Covers. Recommended for individuals without existing skin integrity issues requiring ease of cleaning and heavy usage. Ideal in a Long Term care Setting. 2" cell height.

ADULT SIZES

- IR88LPLTC 15" w x 15" d x 2" h (38cm x 38cm x 5cm)
- IR89LPLTC 15" w x 17" d x 2" h (38cm x 43cm x 5cm)
- IR99LPLTC 17" w x 17" d x 2" h (43cm x 43cm x 5cm)
- IR910LPLTC 17" w x 19" d x 2" h (43cm x 47cm x 5cm)
- IR109LPLTC 19" w x 17" d x 2" h (47cm x 43cm x 5cm)
- IR1010LPLTC 19" w x 19" d x 2" h (47cm x 47cm x 5cm)
- C-ROHO Low-Profile LTC "Custom" (specify size)
- C-ROHO "Custom Design" LTC Cushion (provide diagram)



QUADTRO® SELECT™ -HIGH PROFILE LTC Cushions with Covers.

Recommended for individuals requiring ease of cleaning, heavy usage and positioning capabilities. Ideal in a Long Term Care Setting. 4" cell height.

- QS88LTC 15" w X 15" d X 4" h (38 cm x 38 cm x 10cm)
- QS89LTC 15" w X 17" d X 4" h (38 cm x 43 cm x 10cm)
- QS99LTC 17" w X 17" d X 4" h (43 cm x 43 cm x 10cm)
- QS109LTC 19" w X 17" d X 4" h (47 cm x 43 cm x 10cm)
- QS910LTC 17" w X 19" d X 4" h (43 cm x 47 cm x 10cm)
- QS1010LTC 19" w X 19" d X 4" h (47 cm x 47 cm x 10cm)
- QSC Custom LTC Cushion (please specify size)
- QSC "Custom Design" LTC Cushion (provide diagram)



QUADTRO® SELECT™ -LOW PROFILE LTC Cushions with Covers.

Recommended for individuals without existing skin integrity issues requiring ease of cleaning, heavy usage and positioning capabilities. Ideal in a Long Term Care Setting. 2" cell height.

- QS88LPLTC 15" w X 15" d X 2" h (38 cm x 38 cm x 5 cm)
- QS89LPLTC 15" w X 17" d X 2" h (38 cm x 43 cm x 5 cm)
- QS99LPLTC 17" w X 17" d X 2" h (43 cm x 43 cm x 5 cm)
- QS910LPLTC 17" w X 19" d X 2" h (43 cm x 47 cm x 5 cm)
- QS109LPLTC 19" w X 17" d X 2" h (47 cm x 43 cm x 5 cm)
- QS1010LPLTC 19" w X 19" d X 2" h (47 cm x 47 cm x 5 cm)
- QSC Custom Cushion (please specify size)
- QSC Custom Design LTC Cushion (provide diagram)



ROHO PRO CUSHION.

Combining MERLIN® PROXIMITY SENSOR technology to continuously monitor inflation.

Merlin 1313-Merlin 1815 Merlin 2118

- ROHO PRO SINGLEVALVE for standard high & low profile cushions
- ROHO PRO QUAD SELECT for standard high & low profile cushions
- ROHO PRO SINGLEVALVE CUSTOM for standard high & low profile cushions
- ROHO PRO QUAD SELECT CUSTOM for standard high & low profile cushions



ENHANCER® Cushions (dual valve) with Covers.


Recommended for enhanced pressure distribution, positioning and posture.

ADULT SIZES

- ENH88C 15" w X 15" d (38 cm x 38 cm)
- ENH89C 15" w X 17" d (38 cm x 43 cm)
- ENH99C 17" w X 17" d (43 cm x 43 cm)
- ENH109C 19" w X 17" d (47 cm x 43 cm)
- ENH1010C 19" w X 19" d (47 cm x 47 cm)
- ENH910C 17" w X 19" d (43 cm x 47 cm)
- ENHC Custom Cushion (please specify size)

PEDIATRIC SIZES


- ENH66C 12" w X 12" d (30 cm x 30 cm)
- ENH67C 12" w X 13" d (30 cm x 33 cm)



CONTOUR SELECT® Cushions with Covers.

Recommended for enhanced pressure distribution and posture with the advantage of ISOFLO™

- CS88C 15" w X 15" d (38 cm x 38 cm)
- CS89C 15" w X 17" d (38 cm x 43 cm)
- CS99C 17" w X 17" d (43 cm x 43 cm)
- CS109C 19" w X 17" d (47 cm x 43 cm)
- CS910C 17" w X 19" d (43 cm x 47 cm)
- CS1010 19" w X 19" d (47cm x 47 cm)
- CS Custom Cushion (Please specify size)
- CS Custom Design Cushion (provide diagram)



AIRLITE® Cushions with Removable Cover.

Combines the benefits of sealed air support pad & durable contoured foam base for comfort, enhanced positioning and pressure relief. Designed for client with tactile sensation and who are low risk for pressure ulcers.

- AL 1414 (14" x 14"; 36cm x 36cm)
- AL 1515 (15" X 15"; 38cm X 38cm)
- AL 1616 (16" X 16"; 41cm X 41cm)
- AL 1618 (16" X 18"; 41cm x 46cm)
- AL 1716 (17" x 16"; 43cm x 41cm)
- AL 1717 (17" x 17"; 43 cm x 43 cm)
- AL 1718 (17" x 18"; 43cm x 46cm)
- AL 1720 (17" x 20"; 43cm x 51cm)
- AL 1816 (18" x 16"; 46cm x 41cm)
- AL 1818 (18" x 18"; 46cm x 46cm)
- AL 1820 (18" x 20"; 46cm x 51cm)
- AL 1917 (19" x 17"; 48cm x 43cm)
- AL 1919 (19" x 19"; 48cm x 48cm)
- AL 1920 (19" x 20"; 48cm x 51cm)
- AL 2018 (20" x 18"; 51cm x 46cm)
- ALCUSTMED Airlite Custom Width 14" - 18"/35.5cm-45cm
- ALCUSTLG Airlite Custom Width 19" - 20"/45cm-51cm



NEXUS® SPIRIT™ Cushion with Cover. Combines the stability of a contoured foam base with the therapeutic care for which ROHO® DRY FLOTATION® cushions are known. The contoured foam base is designed to allow for increased stability for transfers and provides exceptional positioning of the pelvis and lower extremities.

- NS1414C 14" w x 14" d (35 cm x 35 cm)
- NS1416C 14" w x 16" d (35 cm x 41 cm)
- NS1515C 15" w x 15" d (38 cm x 38 cm)
- NS1616C 16" w x 16" d (41 cm x 41 cm)
- NS1618C 16" w x 18" d (41 cm x 46 cm)
- NS1717C 17" w x 17" d (43 cm x 43 cm)
- NS1816C 18" w x 16" d (46 cm x 41 cm)
- NS1818C 18" w x 18" d (46 cm x 46 cm)
- NS2018C 20" w x 18" d (51 cm x 46 cm)
- NS CUSTOM DUAL VALVE OR CUSTOM CELL HEIGHT



HARMONY™ Cushion with Cover.

Recommended for individuals seeking moderate pressure relief and enhanced stability.

- HI416C 14" w x 16" d (36 cm X 41 cm)
- HI616C 16" w x 16" d (41 cm x 41 cm)
- HI618C 16" w x 18" d (41 cm x 46 cm)
- HI716C 17" w x 16" d (43 cm x 41 cm)
- HI718C 17" w x 18" d (43 cm x 46 cm)
- HI816C 18" w x 16" d (46 cm x 41 cm)
- HI818C 18" w x 18" d (46 cm x 46 cm)
- H2018C 20" w x 18" d (51 cm x 46 cm)



LTV Cushion with Cover.

The valve allows you to customize the seat to your specific comfort level.

- LTV 18 x 16" Office Wool Cushion
- LTV2B 18 x 16" Ultra Leather - Black




MINI-MAX® Cushions. Recommended for active individuals with low risk of skin breakdown & individuals seeking protection during recreational activities.

- MINI-1515 MINI-MAX* 15" x 15" x 1" (38 cm x 38 cm x 3 cm)
- MINI-1616 MINI-MAX* 16" x 16" x 1" (41 cm x 41 cm x 3 cm)
- MINI-1816 MINI-MAX* 18" x 16" x 1" (46 cm x 41 cm x 3 cm)
- MINI-C "CUSTOM" Mini-Max Cushion



HEAL PAD® Fully adjustable-protects and facilitates healings of heels and elbows

- HEAL PAD Heal pad (single)
- RPL-HPS Replacement heal pad velcro strap
- HP-45 Heal pad pump replacement



REPLACEMENT CELLS/REPAIR KITS

- HP-20 ROHO inflation pump for wheelchair/cushion
- A-5 Rubber patch repair kit
- MRK-I Master repair kit
- CV-1 Inflation valve (with rubber stem)
- G-4 Fastbond 30 NF Adhesive (to replace cells)
- 65010/20 Grommets (Top & Bottom)
- Factory Labour (per repair)

ROHO® HYBRID ELITE™ Cushion with Cover. Combines

the superb skin protection and shape matching ability of ROHO® DRY FLOATATION® with the stability of Jay® base technology. This new cushion design allows for improved pelvic alignment and stability for improved position and function while maintaining optimal immersion and envelopment of both the ischial tuberosities and greater trochanters.




SINGLE VALVE

<input type="checkbox"/>	IRHY1414C	14" w x 14" d (36 cm x 38 cm)
<input type="checkbox"/>	IRHY1415C	14" w x 15" d (36 cm x 38 cm)
<input type="checkbox"/>	IRHY1416C	14" w x 16" d (36 cm x 43 cm)
<input type="checkbox"/>	IRHY1515C	15" w x 15" d (38 cm x 41 cm)
<input type="checkbox"/>	IRHY1516C	15" w x 16" d (38 cm x 41 cm)
<input type="checkbox"/>	IRHY1517C	15" w x 17" d (38 cm x 43 cm)
<input type="checkbox"/>	IRHY1615C	16" w x 15" d (41 cm x 38 cm)
<input type="checkbox"/>	IRHY1616C	16" w x 16" d (41 cm x 41 cm)
<input type="checkbox"/>	IRHY1617C	16" w x 17" d (41 cm x 43 cm)
<input type="checkbox"/>	IRHY1618C	16" w x 18" d (41 cm x 46 cm)
<input type="checkbox"/>	IRHY1716C	17" w x 16" d (43 cm x 41 cm)
<input type="checkbox"/>	IRHY1717C	17" w x 17" d (43 cm x 43 cm)
<input type="checkbox"/>	IRHY1718C	17" w x 18" d (43 cm x 46 cm)
<input type="checkbox"/>	IRHY1816C	18" w x 16" d (46 cm x 41 cm)
<input type="checkbox"/>	IRHY1817C	18" w x 17" d (46 cm x 43 cm)
<input type="checkbox"/>	IRHY1818C	18" w x 18" d (46 cm x 46 cm)
<input type="checkbox"/>	IRHY1820C	18" w x 20" d (46 cm x 51 cm)
<input type="checkbox"/>	IRHY2018C	20" w x 18" d (51 cm x 46 cm)
<input type="checkbox"/>	IRHY2020C	20" w x 20" d (51 cm x 51 cm)
<input type="checkbox"/>	IRHY2218C	22" w x 18" d (56 cm x 46 cm)
<input type="checkbox"/>	IRHY2220C	22" w x 20" d (56 cm x 51 cm)
<input type="checkbox"/>	IRHY2418C	24" w x 18" d (61 cm x 46 cm)
<input type="checkbox"/>	IRHY2420C	24" w x 20" d (61 cm x 51 cm)

DUAL VALVE

<input type="checkbox"/>	2RHY1414C	14" w x 14" d (36 cm x 38 cm)
<input type="checkbox"/>	2RHY1415C	14" w x 15" d (36 cm x 38 cm)
<input type="checkbox"/>	2RHY1416C	14" w x 16" d (36 cm x 43 cm)
<input type="checkbox"/>	2RHY1515C	15" w x 15" d (38 cm x 41 cm)
<input type="checkbox"/>	2RHY1516C	15" w x 16" d (38 cm x 41 cm)
<input type="checkbox"/>	2RHY1517C	15" w x 17" d (38 cm x 43 cm)
<input type="checkbox"/>	2RHY1615C	16" w x 15" d (41 cm x 38 cm)
<input type="checkbox"/>	2RHY1616C	16" w x 16" d (41 cm x 41 cm)
<input type="checkbox"/>	2RHY1617C	16" w x 17" d (41 cm x 43 cm)
<input type="checkbox"/>	2RHY1618C	16" w x 18" d (41 cm x 46 cm)
<input type="checkbox"/>	2RHY1716C	17" w x 16" d (43 cm x 41 cm)
<input type="checkbox"/>	2RHY1717C	17" w x 17" d (43 cm x 43 cm)
<input type="checkbox"/>	2RHY1718C	17" w x 18" d (43 cm x 46 cm)
<input type="checkbox"/>	2RHY1816C	18" w x 16" d (46 cm x 41 cm)
<input type="checkbox"/>	2RHY1817C	18" w x 17" d (46 cm x 43 cm)
<input type="checkbox"/>	2RHY1818C	18" w x 18" d (46 cm x 46 cm)
<input type="checkbox"/>	2RHY1820C	18" w x 20" d (46 cm x 51 cm)
<input type="checkbox"/>	2RHY2018C	20" w x 18" d (51 cm x 46 cm)
<input type="checkbox"/>	2RHY2020C	20" w x 20" d (51 cm x 51 cm)
<input type="checkbox"/>	2RHY2218C	22" w x 18" d (56 cm x 46 cm)
<input type="checkbox"/>	2RHY2220C	22" w x 20" d (56 cm x 51 cm)
<input type="checkbox"/>	2RHY2418C	24" w x 18" d (61 cm x 46 cm)
<input type="checkbox"/>	2RHY2420C	24" w x 20" d (61 cm x 51 cm)

SPECIAL APPLICATION CUSHIONS. Unique design solutions for every seating problem.




ADAPTOR® 9" x 13" x 1"
sealed 1 square inch air cell
padding*
*can be cut to fit virtually
any surface
* DOES NOT INCLUDE COVER



* DOES NOT INCLUDE COVER

- SACRAL82 Sacral Support Cushion (4"h. x 15"w) (10 cm x 38 cm)
- SACRAL92 Sacral Support Cushion (4"h. x 17"w) (10 cm x 42 cm)
- SACRAL102 Sacral Support Cushion (4"h. x 19"w) (10 cm x 47 cm)



* DOES NOT INCLUDE COVER

Shower CommodesCushion*
(Low profile - 2" cell height only)

- SEAT CMD 88 (15" x 15") (38 cm x 38 cm)
- SEAT CMD 89 (15" x 17") (38 cm x 43 cm)
- SEAT CMD 99 (17" x 17") (43 cm x 43 cm)
- SEAT CMD 910 (17" x 18") (43 cm x 47 cm)
- SEAT CMD CUSTOM

RECLINER-C
Recliner wheelchair cushion*
(Made to fit per order, up to 230 cells overall)

CUST RECLINER-C "CUSTOM"
Recliner wheelchair cushion*
PRICE UPON REQUEST.



T-SEAT Toilet seat cushion*
Fits standard toilet seat
(no custom size available)
* DOES NOT INCLUDE COVER



PACK -IT®-B
with cover, design to suit short
term needs such as bathing or
travelling 16" x 9" x 1"
41 cm x 23 cm x 3 cm

LUM-C Lumbar support*
cushion* (6.5"h x 15"w)
(17 cm x 38 cm)

* DOES NOT INCLUDE COVER

Specify attachment:
 Snap Velcro Strap
Specify cell height: 1" to 4"



LATERAL
Lateral support cushion* (please specify).....CUSTOM
* Does not include cover

SEATING ACCESSORIES AND COVERS




MERLIN® PROXIMITY SENSOR. Precisely and continually monitors your ROHO Cushion's inflation level, keeping you worry free.

- Merlin1313 Fits cushion size 7 cells x 7 cells to 8 cells x 10 cells
- Merlin1615 Fits cushion size 9 cells x 8 cells to 9 cells x 11 cells
- Merlin1815 Fits cushion size 10 cells x 8 cells to 12 cells x 12 cells
- Merlin2118 Fits cushion size 12 cells x 10 cells to 12 cells x 12 cells



CUSHION RETAINER Keep cushion from sliding out of the back of the chair.


- CR-1 Cushion Retainer



SOLID SEAT INSERTS



• Custom sizes available
Not recommended for
individuals over 340 lbs.



CONTOUR BASE

(used with sling seat)

- BASE1515 15" X 15" (38 cm x 38 cm)
- BASE1715 17" x 15" (43 cm x 38 cm)
- BASE1717 17" x 18" (43 cm x 45 cm)

PLANAR SOLID SEAT INSERT

(used with sling seat)

- SI 1416 14"w x 16"d (36cm x 41 cm)
- SI 1616 16"w x 16"d (41 cm x 41 cm)
- SI 1618 16"w x 18"d (41 cm x 46 cm)
- SI 1818 18"w x 18"d (46cm x 46 cm)
- SI 1820 18"w x 20"d (46 cm x 51 cm)
- SI 2020 20"w x 20"d (51 cm x 51 cm)
- SI Custom (please specify)

NON-SLIP Replacement Cushion Covers

Assist in slide transfers and minimizes wear on the cushion.

Universal Cover Standard sizes (for single, dual valve, Quadro Select, Enhancer & Contour Select)

STANDARD SIZES (Match with cushion size listed below) **CUSTOM CUSHIONS** (Any other than those listed)

NEXUS SPIRIT STD SIZES

HARMONY STD SIZES

AIRLITE STD SIZES

LTV

- SCLTVW - WOOL
- SCLTVL - BLACK ULTRA LEATHER

HIGH or MID PROFILE

- UC88
- UC89
- UC99
- UC910
- UC109
- UC1010

LOW PROFILE

- UC88LP
- UC89LP
- UC99LP
- UC910LP
- UC109LP
- UC1010LP

MINI MAX

- SCMT15
- SCMT16
- SCMT18

NEXUS

- SCNS1414
- SCNS1416
- SCNS1515
- SCNS1616
- SCNS1618
- SCNS1717
- SCNS1816
- SCNS1818
- SCNS2018

HARMONY

- SCH1416
- SCH1616
- SCH1618
- SCH1716
- SCH1718
- SCH1816
- SCH1818
- SCH2018

AIRLITE

- | | |
|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> SCAL1414 | <input type="checkbox"/> SCAL1816 |
| <input type="checkbox"/> SCAL1515 | <input type="checkbox"/> SCAL1818 |
| <input type="checkbox"/> SCAL1616 | <input type="checkbox"/> SCAL1820 |
| <input type="checkbox"/> SCAL1618 | <input type="checkbox"/> SCAL1917 |
| <input type="checkbox"/> SCAL1716 | <input type="checkbox"/> SCAL1919 |
| <input type="checkbox"/> SCAL1717 | <input type="checkbox"/> SCAL1920 |
| <input type="checkbox"/> SCAL1718 | <input type="checkbox"/> SCAL2018 |
| <input type="checkbox"/> SCAL1720 | |

LTC CUSHION COVERS

Protect your ROHO cushion from incontinent incidents. Cover material protects against microfungi and bacteria. Allows for easy cleaning to ensure a clean, sterile sitting environment for the user. (for single, dual valve, Quadro Select, Enhancer & Contour Select)

- | | | | | |
|--------------------------------|---------------------------------|----------------------------------|----------------------------------|------------------------------------|
| <input type="checkbox"/> LTC66 | <input type="checkbox"/> LTC98 | <input type="checkbox"/> LTC911 | <input type="checkbox"/> LTC66LP | <input type="checkbox"/> LTC98LP |
| <input type="checkbox"/> LTC67 | <input type="checkbox"/> LTC99 | <input type="checkbox"/> LTC1010 | <input type="checkbox"/> LTC67LP | <input type="checkbox"/> LTC99LP |
| <input type="checkbox"/> LTC77 | <input type="checkbox"/> LTC108 | <input type="checkbox"/> LTC1011 | <input type="checkbox"/> LTC77LP | <input type="checkbox"/> LTC108LP |
| <input type="checkbox"/> LTC78 | <input type="checkbox"/> LTC109 | <input type="checkbox"/> LTC1110 | <input type="checkbox"/> LTC78LP | <input type="checkbox"/> LTC109LP |
| <input type="checkbox"/> LTC87 | <input type="checkbox"/> LTC119 | | <input type="checkbox"/> LTC87LP | <input type="checkbox"/> LTC810LP |
| <input type="checkbox"/> LTC88 | <input type="checkbox"/> LTC810 | | <input type="checkbox"/> LTC88LP | <input type="checkbox"/> LTC910LP |
| <input type="checkbox"/> LTC89 | <input type="checkbox"/> LTC910 | | <input type="checkbox"/> LTC89LP | <input type="checkbox"/> LTC1010LP |



HYBRID ELITE™ CUSHION COVERS.

Two-way stretch cover with encased comfort foam cells, conforms to the individual's body shape. Solid, non-skid base material offers a durable easy to clean surface.

Standard sizes

- | | |
|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> UCHY1414 | <input type="checkbox"/> UCHY1816 |
| <input type="checkbox"/> UCHY1415 | <input type="checkbox"/> UCHY1817 |
| <input type="checkbox"/> UCHY1416 | <input type="checkbox"/> UCHY1820 |
| <input type="checkbox"/> UCHY1515 | <input type="checkbox"/> UCHY2018 |
| <input type="checkbox"/> UCHY1516 | <input type="checkbox"/> UCHY2020 |
| <input type="checkbox"/> UCHY1517 | <input type="checkbox"/> UCHY2218 |
| <input type="checkbox"/> UCHY1615 | <input type="checkbox"/> UCHY2220 |
| <input type="checkbox"/> UCHY1616 | <input type="checkbox"/> UCHY2418 |
| <input type="checkbox"/> UCHY1617 | <input type="checkbox"/> UCHY2420 |
| <input type="checkbox"/> UCHY1618 | |
| <input type="checkbox"/> UCHY1717 | |
| <input type="checkbox"/> UCHY1718 | |

HYBRID ELITE™ HEAVY DUTY CUSHION COVERS.

Protect your Hybrid Elite Roho cushion from incontinent incidents. Cover material protects against microfungi and bacteria. Allow for easy cleaning to ensure a clean, sterile sitting environment (for single & dual valve Hybrid Elite cushions). Standard sizes

- | | |
|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> HDHY1414 | <input type="checkbox"/> HDHY1816 |
| <input type="checkbox"/> HDHY1415 | <input type="checkbox"/> HDHY1817 |
| <input type="checkbox"/> HDHY1416 | <input type="checkbox"/> HDHY1820 |
| <input type="checkbox"/> HDHY1515 | <input type="checkbox"/> HDHY2018 |
| <input type="checkbox"/> HDHY1516 | <input type="checkbox"/> HDHY2020 |
| <input type="checkbox"/> HDHY1517 | <input type="checkbox"/> HDHY2218 |
| <input type="checkbox"/> HDHY1615 | <input type="checkbox"/> HDHY2220 |
| <input type="checkbox"/> HDHY1616 | <input type="checkbox"/> HDHY2418 |
| <input type="checkbox"/> HDHY1617 | <input type="checkbox"/> HDHY2420 |
| <input type="checkbox"/> HDHY1618 | |
| <input type="checkbox"/> HDHY1717 | |
| <input type="checkbox"/> HDHY1718 | |

NOTE 1: Ventilation/drainage openings through base of cushion available. See Custom Price
 NOTE 2: All cushion prices include a cover, inflation pump, repair kit, and operating instructions, except where noted by an *
 NOTE 3: All measurements for cushions are approximate from edge to edge.

SIZING CHART - ROHO® DRY FLOATATION® Cushions

The following sizes are consistent with all ROHO DRY FLOATATION standard seat cushions, excluding the NEXUS® SPIRIT™ and AirLITE® cushion. The number of cells are consistent with model numbers listed next to each product. Wheelchair seat size should be used to determine the appropriate number of cells required to fit your seat.

For the NEXUS® SPIRIT™ and AirLITE cushions, the model number indicates the dimensions of the wheelchair seat in inches, width x length.

Most ROHO cushions can be customized beyond standard specifications. For pricing on custom cushions, call 800-265-9830
 Custom Sizes greater than 160 cells - Priced Upon Request

TO FIT CHAIR SIZE LISTED BELOW	SELECT THIS NUMBER OF CELLS
33 cm x 33 cm / 13" x 13"	7 wide x 7 deep
36 cm x 36 cm / 14" x 14"	8 wide x 8 deep*
33 cm x 38 cm / 13" x 15"	7 wide x 8 deep
38 cm x 33 cm / 15" x 13"	8 wide x 7 deep
38 cm x 38 cm / 15" x 15"	8 wide x 8 deep
38 cm x 41 cm / 15" x 16"	8 wide x 9 deep
38 cm x 48 cm / 15" x 18"	8 wide x 10 deep
41 cm x 38 cm / 16" x 15"	9 wide x 8 deep*
41 cm x 41 cm / 16" x 16"	9 wide x 9 deep*
43 cm x 43 cm / 17" x 17"	9 wide x 9 deep
41 cm x 48 cm / 16" x 18"	9 wide x 10 deep*
48 cm x 38 cm / 18" x 15"	10 wide x 8 deep*
48 cm x 41 cm / 18" x 16"	10 wide x 9 deep*
48 cm x 48 cm / 18" x 18"	10 wide x 10 deep*

THERAPEUTIC SURFACES

Roho Mattress Sections

- 3400 Roho Neoprene Section (20" x 34" / 51cm x 86cm)
 - 3800-1 Prodigy Mattress Overlay Single Section (27" x 35" / 68cm x 89cm)
 - 3800-3 Prodigy Full Mattress Overlay (80" x 35" / 203cm x 89cm)
 - UC MATT Cover for Neoprene section
 - IC MATT Incontinent Cover for Neoprene Section
- Custom Neoprene Section available (Call for Quote)

Leveling Pads

- RLP3480 Roho Leveling Pad 34" x 80" (with pre-cut centre) 86cm x 203cm
- RLP2134 Roho Leveling Pad 21" x 34" / 54cm x 86cm
- RLP3134 Roho Leveling Pad 31" x 34" / 79cm x 86cm
- PLP3580 Prodigy Leveling Pad 35" x 80" / 89cm x 203cm (with pre-cut centre)
- PLP2735 Prodigy Leveling Pad 27" x 35" / 68cm x 89cm
- PLP2722 Prodigy Leveling Pad 27" x 22" / 68cm x 56cm

Mattress Covers

- RRCOVER3 3" Deep Reusable Cover for Roho Sections 36" x 81"/91cm x 206cm
- RRCOVER9 9" Deep Reusable Cover for Roho Sections and or Mattress
- 988 Maxiflo Underpad 23" x 36" / 58cm x 91cm pkg. 12

Vitacare Medical Products V4 Sectional Mattresses

The Vitacare V4 Mattress system with Roho creates the optimal therapeutic environment for skin healing and protection.



Roho Neoprene Insert is 19" Wide

- 3040NU V4 with 19" Opening
- 3040NU-219 V4 with 2 19" Openings
- 3040NU-138 V4 with 1 38" Opening
- 3040NU-146 V4 with 2 19" Openings & moveable divider
- 3040NU-157 V4 with 1 57" Opening
- 304I 19" Foam insert to fit into 3040NU
- 306I 27" Foam Insert to fit into 3060NU
- 3004 V4 Recovery 5 Sectional Cover specify openings
- Custom Mattress Call Vitacare for Quote 1-800-263-9068

- 3080NU V8 Mattress 36" x 80" x 6' / 91cm x 203cm x 15 cm
Non-powered Mattress replacement designed for individuals who are at risk for breakdown and/or have a Stage I or Stage II ischemic ulcer.
Zoned for maximum immersion in all areas of the body.



SelectAir™ Low Air Loss Therapy Systems

- SASTDSYSK001 **SelectAir MATTRESS SYSTEM STANDARD**
System Includes: *SelectAir* pump (SASTDPUMP), mattress with foot pillow (SASTDMATT), zippered cover with foot pillow (SASTDCOV) and a 1" foam pad (SASTDFOAMI).
- SAMAXSYSK004 **SelectAir MAX MATTRESS SYSTEM**
System Includes: *SelectAir* Max pump (SAMAXPUMP), mattress with foot pillow (SASTDMATT), zippered cover with foot pillow (SASTDCOV) and a 1" foam pad (SASTDFOAMI).
- SAMAXSYSK010 **SelectAir MAX NON-FOOTED MATTRESS SYSTEM**
System Includes: *SelectAir* Max pump (SAMAXPUMP), non-footed mattress (SANOFTMATT), non-footed zippered cover (SANOFTCOV), and a 1" foam pad (SASTDFOAMI).
- SASTDSYSK008 **SelectAir NON-FOOTED MATTRESS SYSTEM**
System Includes: *SelectAir* pump (SASTDPUMP), non-footed mattress (SANOFTMATT), non-footed zippered cover (SANOFTCOV), and a 1" foam pad (SASTDFOAMI).
- SAOVERLAYSYS **SelectAir OVERLAY SYSTEM STANDARD**
System Includes: *SelectAir* pump (SASTDPUMP), 4" overlay mattress (SAOVERLAYMATT) and overlay mattress sheet (SAOVERLAYSHEET).
- SABARISYS39 **BariSelect 39" BARIATRIC MATTRESS SYSTEM**
System Includes: *BariSelect* pump (SABARIPUMP), 39" mattress (SA39MATT), 39" zippered cover (SA39COV) and a 35" foam pad (SABARIFOAMI).
- SABARISYS42 **BariSelect 42" BARIATRIC MATTRESS SYSTEM**
System Includes: *BariSelect* pump (SABARIPUMP), 42" mattress (SA42MATT), 42" zippered cover (SA42COV) and a 35" foam pad (SABARIFOAMI).
- SABARISYS48 **BariSelect 48" BARIATRIC MATTRESS SYSTEM**
System Includes: *BariSelect* pump (SABARIPUMP), 48" mattress (SA48MATT), 48" zippered cover (SA48COV) and a 35" foam pad (SABARIFOAMI).
- SABARISYS54 **BariSelect 54" BARIATRIC MATTRESS SYSTEM**
System Includes: *BariSelect* pump (SABARIPUMP), 54" mattress (SA54MATT), 54" zippered cover (SA54COV) and a 35" foam pad (SABARIFOAMI).
- SABARISYS60 **BariSelect 60" BARIATRIC MATTRESS SYSTEM**
System Includes: *BariSelect* pump (SABARIPUMP), 60" mattress (SA60MATT), 60" zippered cover (SA60COV) and a 35" foam pad (SABARIFOAMI).
- SAHYOVMATT **Hybrid Select™ LOW AIR LOSS SYSTEM**
System Includes: Mattress overlay, pump, cover
Heal Float Cushion with cover, instruction booklet

SELECT PROTECT™ SAFETY ENHANCEMENT SYSTEM

SelectProtect is designed to be used exclusively with standard ROHO SelectAir or SelectAir Max Low Air Loss Mattress Replacement Systems to increase patient bed safety. Integrated side bolsters help provide a safe sleeping environment for patients. SelectProtect is designed to be used without bed side rails*. It is ideal for use in many healthcare settings where bed side rails are not used including long-term care facilities and in the home.



- | | | |
|--------------------------|----------------|--|
| <input type="checkbox"/> | SASPSTDSYSCO I | Select Air Matt System Standard w/ SA Protect Std. Cover |
| <input type="checkbox"/> | SASPMAXSYSCO I | Select Air Max Matt System w/ SA Protect Std Cover |
| <input type="checkbox"/> | SASPMAXSYSCO2 | Select Air Max Matt System w/ SA Protect no foot cover |
| <input type="checkbox"/> | SAPROTECTSTD | Cover, Select Protect std complete |
| <input type="checkbox"/> | SAPROTECTNOFT | Cover, Select Protect no foot complete |
| <input type="checkbox"/> | SASTDCOVSP | Cover, "ONLY" Select Protect Std cover |
| <input type="checkbox"/> | SANOFTCOVSP | Cover, "ONLY" Select Protect no foot |
| <input type="checkbox"/> | SAPROTSHBOL | Replacement shoulder bolster (foam) |
| <input type="checkbox"/> | SAPROTFTBOL | Replacement foot bolster (foam) |

PLEASE NOTE: * Use of SelectProtect should be determined based upon facility protocols, applicable laws, regulations, guidelines and individual user safety considerations. It does not replace the need for patient fall prevention measures.



SelectAir™ Therapy System Replacement Parts

- | | | |
|--------------------------|---------------|---------------------------------------|
| <input type="checkbox"/> | SASTDPUMP | SelectAir pump |
| <input type="checkbox"/> | SAMAXPUMP | SelectAir Max pump |
| <input type="checkbox"/> | SABARIPUMP | SelectAir BariSelect pump |
| <input type="checkbox"/> | SAHYBMATT | SelectAir Hybrid Mattress |
| <input type="checkbox"/> | SASTDMATT | Standard mattress with foot pillow |
| <input type="checkbox"/> | SASHORTMATT | Short mattress with foot pillow |
| <input type="checkbox"/> | SA2SHORTMATT | Short/Short mattress with foot pillow |
| <input type="checkbox"/> | SANOFTMATT | Standard non-footed mattress |
| <input type="checkbox"/> | SAOVERLAYMATT | Overlay mattress |
| <input type="checkbox"/> | SA39MATT | 39" mattress (84" x 39" x 8") |
| <input type="checkbox"/> | SA42MATT | 42" mattress (84" x 42" x 8") |
| <input type="checkbox"/> | SA48MATT | 48" mattress (84" x 48" x 8") |
| <input type="checkbox"/> | SA54MATT | 54" mattress (84" x 54" x 8") |
| <input type="checkbox"/> | SA60MATT | 60" mattress (84" x 60" x 8") |

continued

<input type="checkbox"/>	SASTDCOV	Standard cover with foot pillow
<input type="checkbox"/>	SASHORTCOV	Short cover with foot pillow
<input type="checkbox"/>	SANOFTCOV	Standard non-footed cover
<input type="checkbox"/>	SA2PCCOV	2 piece cover with foot pillow
<input type="checkbox"/>	SA30COV	39" mattress cover
<input type="checkbox"/>	SA42COV	42" mattress cover
<input type="checkbox"/>	SA48COV	48" mattress cover
<input type="checkbox"/>	SA54COV	54" mattress cover
<input type="checkbox"/>	SA60COV	60" mattress cover
<input type="checkbox"/>	SA2PCCOVBTM	Bottom for 2 piece cover with foot pillow
<input type="checkbox"/>	SA2PCCOVTOP	Top for 2 piece cover with foot pillow
<input type="checkbox"/>	SAMATTSHEET	Sheet for mattress with foot pillow
<input type="checkbox"/>	SAOVERLAY SHEET	Overlay mattress sheet
<input type="checkbox"/>	HYFOAM BASE	2" Hybrid Foam Base
<input type="checkbox"/>	SAHYBSHEET	Select Air Hybrid Sheet
<input type="checkbox"/>	SASTDFOAMI	1" foam pad
<input type="checkbox"/>	SABARIFOAM	Barifoam pad 35"
<input type="checkbox"/>	SABAG	System carry bag
<input type="checkbox"/>	SABARIBAG	Bariatric system carry bag
<input type="checkbox"/>	SACONTROL	Hand control for deluxe/Bari
<input type="checkbox"/>	SAFILTER	Filter, package of 6
<input type="checkbox"/>	SAHANDLE	Handle
<input type="checkbox"/>	SAHOSE	Hose
<input type="checkbox"/>	SAPOWERCORD	Power cord

THERAPEUTIC PILLOWS

Sleepmatterzz®

Sleepmatterzz Original Pillow - Box Cover (Batting on one side)

<input type="checkbox"/>	SMPIRCRB	Single Valve Standard 14" x 19"/36cm x 48cm
<input type="checkbox"/>	SMP2RCRB	Dual Valve Standard 14" x 19"/36cm x 48cm
<input type="checkbox"/>	SMPIRCRBQ	Single Valve Queen 14" x 25"/36cm x 64cm
<input type="checkbox"/>	SMP2RCRBQ	Dual Valve Queen 14" x 25"/36cm x 64cm

Sleepmatterzz Single Valve Classic Pillow - Classic Cover (Batting on two sides)

<input type="checkbox"/>	SMPIRCRCLS	Single Valve Classic Standard 14"x19"/36cmx48cm
<input type="checkbox"/>	SMP2RCRCLS	Dual Valve Classic Standard 14"x19"/36cmx48cm
<input type="checkbox"/>	SMPIRCRCLSQ	Single Valve Classic Queen 14"x25"/36cmx64cm
<input type="checkbox"/>	SMP2RCRCLSQ	Dual Valve Classic Queen 14"x25"/36cmx64cm
<input type="checkbox"/>	SMPIRCRCLSK	Single Valve Classic King 14"x29"/36cmx74cm
<input type="checkbox"/>	SMP2RCRCLSK	Dual Valve Classic King 14"x29"/36cmx74cm

TERMS AND POLICIES

1. RETURN POLICY: ALL RETURNS REQUIRE PRIOR AUTHORIZATION AND ARE SUBJECT TO A RESTOCKING CHARGE.
2. REPAIR CHARGES: A REPAIR CHARGE OF \$75.00 PLUS SHIPPING WILL BE CHARGED FOR ALL CUSHIONS AND NON-POWERED MATTRESS PRODUCTS NOT COVERED UNDER WARRANTY.
WE RESERVE THE SOLE RIGHT TO DETERMINE WHETHER CUSHION OR MATTRESS IS REPAIRABLE. COMPLEX REPAIRS MAY REQUIRE ADDITIONAL CHARGE. POWER SURFACES CALL FOR QUOTE.
3. ALL CUSHION AND MATTRESS PRICES INCLUDE CUSHION, MATTRESS, INFLATION PUMP, RUBBER PATCH, REPAIR KIT, AND OPERATING INSTRUCTIONS. CUSHION PRICES ALSO INCLUDE COVER UNLESS NOTED BY *
4. FREIGHT: ALL MATTRESS AND PILLOW PRODUCTS F.O.B. LACHINE, QUEBEC.

ADVANCED HEALTH CARE PRODUCTS reserve the right to make changes to prices and specifications without notice.

Other Products Available From Advanced Health Care

 **Bodypoint**[®]

XSENSOR[®] Technology Corporation
Innovators in Pressure Imaging



showerbuddy[™]



PRIZM MEDICAL INC.

Carex[™]
caring for you

ACCESS[™]

VITAL

Thomashilfen

DISTRIBUTED IN CANADA BY:

AHC
ADVANCED HEALTH CARE PRODUCTS G.P.
Enhancing Your Quality of Life

2221 - 46th Ave., Lachine, Quebec H8T 3C9
Tel: 514-636-7575 • 1-800-265-9830
Fax: 877-636-8944
www.advancedhealthcare.ca